



TCS WASHING FASTNESS SAMPLE DYEING MACHINE



Truly we fulfill your wish

LABORATORY WASHING FASTNESS:

Laboratory Washing fastness tester is an ideal machine for shade matching and for testing of various types of dyes up to the dyeing temperature of 950c. In this machine both substrate and dyeing liquor is continuously moving resulting in uniform dyeing results.

Sample for dyeing and ready dyeing solution are put in s.s beaker and all prepared beakers are inserted in to the carriage. Carriage rotates at 24 rpm in a heated water bath. A micro processor controller PC-07 is used for time, temperature control. Machine is available in following models.

08 x 500ml, 12 x 500ml,

Technical Specifications:

Carriage rpm	: 40 rpm.
CARRIAGE ROTATION	: FORWARD.
Automation	: through PC-07.
Electrical Heating for Temperature Up	: To 98°C
Heating	: 3kw.
Dimensions Approx	: L 750mmxW330mmxH450mm,
Electrical Supply	: 3kw , 215V single phase, 50Hz

Factory : Global Industrial Park, A-14, Survey No: 26/2, Near Nahuli Railway Crossing
Off. N.H. 8, Valwada, VAPI – 396195, Dist: Valsad, Gujarat. E-mail:
tcsindustries2012@yahoo.com, Mob: 09998642044, 09173480650.

RESIDENCE: Plot No: 395, Neela Co-Op.Hsg.Society.Ltd, Near Ambamata Mandir,
Koparli Road, Gunjan, G.I.D.C, VAPI-396195. Mob: 09998642044, 09173480650.